

BORAZ

## NEW BORA 2 Iconic reliability

The Bora turbine has continued to reaffirm the quality of its engineering and its high level of reliability, both of which have remained exemplary. It is now Bien-Air's bestselling turbine. Because of the trust you have placed in it, we wanted to ensure that this legacy continues

SWISS 🛨 MADE

### IN 1991 an icon was born

#### From now on, the past speaks for itself; the future is guaranteed.

BORAZ

The reliability and service life of the Bora turbine are self-evident.



The new design of the Bora 2 turbine provides even better performance, whilst retaining the legendary reliability which has characterised previous versions since 1991.

SWISS 📑 MADE

·C No





#### **EXEMPLARY RELIABILITY** AND ENHANCED PERFORMANCE.

The new design of the Bora 2 offers even better performance, whilst retaining the legendary reliability which has characterised previous versions since 1991.

#### MORE ROBUST.

The head of the Bora 2 is made entirely from stainless steel. This increases its durability, making it more resistant to impacts.

#### MORE **POWERFUL.**

With a **20 %** increase in both power and torque, the Bora 2 allows you to cut **twice as fast.** Save precious time!

# **NEW DESIGN.**



SWISS 🕂 MADE

#### NOTHING CAN HINDER YOU.

To ensure better irrigation of the area being treated, the Accu-Spray Quattro Mix<sup>™</sup> system has four asymmetrical air/water sprays which converge on the end of the bur, with laser precision. In combination with the LED lighting, perfect visibility of the oral cavity is guaranteed.











#### EXTENDED SERVICE LIFE.

The optimised design of the rotor and the new ceramic ball bearings with enhanced wear-resistance guarantee an extended service life **35 %** longer than the previous generation.

#### INCREASED PRECISION.

The Accu-Chuck PreciPlus chuck design ensures the most secure clamping of the bur. Noise is significantly reduced, thereby lending the Bora 2 unequalled practitioner and patient comfort.

#### MORE SAFETY.

Integrated into the push-button in the handpiece head, the patented **Cool Touch** heat-guard technology prevents a build-up of heat and reduces the risk of burns to your patients.



# **SMALLER IS BETTER.**





#### REDUCED RISKS OF CROSS CONTAMINATION.

The only handpiece with an sterilizable **anti-retraction valve** is combined with the only coupling containing exhaust flow anti-retraction.

Studies have proven anti-retraction valves considerably reduce the flow of fluids from the oral cavity back into the dental unit's hoses. Alongside this, the protective **Sealed Head** mechanism reduces the aspiration of fluids into the instrument head.

#### SMALLER HEAD.

Thanks to its smaller head, which is less than **12 mm** in diameter and just **13 mm** long, the Bora 2 increases your field of vision and allows better access to the back of the mouth, for enhanced comfort.



#### BORA 2 LK

20 watts 59 dBA Head size Ø 11.9 mm 13.2 mm 4 asymmetrical mixed sprays anti-retraction valve LED Light 2 years warranty

Bora 2 LK REF 1601153-001



#### WWW.BIENAIR.COM

SWISS 🕂 MADE



Read more